Comparative study between two automated urine strip analyzers: Miditron M and Uriscan Pro+

Prapawadee Ekawong*
Viroj Wiwanitkit**

Ekawong P, Wiwanitkit V. Comparative study between two automated urine strip analyzers: Miditron M and Uriscan Pro+. Chula Med J 2007 Sep - Oct; 51(9): 405 - 10

Background

Presently, due to the advancement of technology in laboratory medicine, many automated analyzers have been launched for urinalysis. Urine strip analyzer or urine strip reader is an automated urine analyzer for semi-quantitative analysis of urine biochemistry. The machine has been widely used in health check up. Here, the authors reported a comparative study between two automated urine strip analyzer machines from two different companies used at King Chulalongkorn Memorial Hospital.

Materials and Methods

This study was comparative between the two automated urine strip analyzers: Miditron M and Uriscan Pro+. One hundred and ninety-nine individual urine samples were included. For each subject, urine strip tests by both methods were performed.

Results

The agreement of glucose, bilirubin, pH and specific gravity determined by both machines were: 100 %, 100 %, 70.8 % and 65.3 %, respectively. A significant correlation was found when Miditron M results of pH and specific gravity were compared with those obtained from Uriscan Pro⁺.

^{*} Division of Laboratory Medicine, King Chulalongkorn Memorial Hospital

^{**} Department of Laboratory Medicine, Faculty of Medicine, Chulalongkorn University

Conclusion

According to this study, the author found that although

the two machines did provide 100 % agreement of glucose

and bilirubin results from both machines, and the agreement

of pH and specific gravity determined by both machines

equals 70.8 % and 65.3 %, respectively, a significant

correlation between the two machines was found. The authors

hereby concluded that using either urine strip analyzer

machine can provide similar analytical results of glucose,

bilirubin, pH and specific gravity.

Keywords

Urinalysis, Automated, Miditron M, Uriscan Pro+

Reprint request: Ekawong P. Division of Laboratory Medicine, King Chulalongkorn Memorial

Hospital, Bangkok 10330, Thailand.

Received for publication. March 15,2005.

ประภาวดี เอกวงศ์, วิโรจน์ ไววานิชกิจ. การศึกษาเปรียบเทียบเครื่องตรวจแถบตรวจปัสสาวะ Miditron M และ Uriscan Pro+. จุฬาลงกรณ์เวชสาร 2550 ก.ย. - ต.ค; 51(9): 405 - 10

พื้นฐาน

ปัจจุบันมีความก้าวหน้าเกี่ยวกับเทคโนโลยีเครื่องอัตโนมัติอย่างมากใน ด้านอายุรกรรมห้องปฏิบัติการ มีการนำเครื่องตรวจวิเคราะห์มาใช้ในการ ตรวจปัสสาวะ เครื่องตรวจแถบตรวจปัสสาวะเป็นเครื่องอัตโนมัติชนิดหนึ่ง สำหรับตรวจวิเคราะห์ปัสสาวะแบบคุณภาพวิเคราะห์ด้านชีวเคมี ได้มีการ นำเครื่องดังกล่าวมาใช้อย่างกว้างขวางในงานด้านการตรวจสุขภาพ ใน การศึกษานี้ผู้นิพนธ์ได้ทำการศึกษาเปรียบเทียบเครื่องตรวจแถบตรวจ ปัสสาวะ 2 ชนิดคือ Miditron M และ Uriscan Pro+.

วัสดุและวิธีการ

ได้ทำการเปรียบเทียบการวิเคราะห์ด้วยเครื่อง Miditron M และ Uriscan Pro+ โดยใช้ตัวอย่างปัสสาวะจากตัวอย่างส่งตรวจ ที่ส่งตรวจยังหน่วยตรวจ ปัสสาวะจำนวน 199 ตัวอย่าง

ผลการศึกษา

ค่าความสอดคล้องของผลการตรวจวิเคราะห์สำหรับ glucose, bilirubin, pH และ specific gravity จากเครื่องทั้งสองเท่ากับ 100 %, 100 %, 70.8 % และ 65.3 % ตามลำดับ และพบว่า ค่า pH และ ค่า specific gravity ที่ตรวจวิเคราะห์จากเครื่องตรวจแถบตรวจปัสสาวะ Miditron M มีความ สัมพันธ์กับเครื่อง Uriscan Pro+

สรุป

พบว่าเครื่องทั้งสองให้ผลการตรวจที่สอดคล้องกันสูง 100 % สำหรับการ
วิเคราะห์ glucose และ bilirubin ดังนั้นการใช้เครื่องตรวจชนิดใดก็ตาม
ย่อมให้ผลการตรวจที่เทียบเคียงกันได้ในกรณีการวิเคราะห์ glucose และ
bilirubin และแม้ว่าค่าความสอดคล้องของผลการตรวจวิเคราะห์สำหรับค่า
pH และ ค่า specific gravity จะเท่ากับ 70.8 % และ 65.3 % ตามลำดับ
แต่ก็มีความสัมพันธ์กัน ดังนั้นการใช้เครื่องตรวจวิเคราะห์แถบทดสอบ
ปัสสาวะชนิดใดก็ตาม ย่อมให้ผลการตรวจที่เทียบเคียงกันได้ในกรณีการ
ตรวจวิเคราะห์ glucose bilirubin pH และ specific gravity

คำสำคัญ

การตรวจปัสสาวะ, อัตโนมัติ, Miditron M, Uriscan Pro+

Urinary analysis (UA) is a basic laboratory investigation in laboratory medicine. It is useful for screening and detection of many diseases. (1,2) In addition, it is recommended for a routine annual check up test. (3) Generally, a routine UA-is composed of two parts: microscopic urine sedimentation analysis and urine biochemistry test by urine strip. (1,2)

Presently, due to the advancement of technology in laboratory medicine, many automated analyzers have been launched for urinalysis. (4) Urine strip analyzer or urine strip reader is an automated urine analyzer for semi-quantitative analysis of urine biochemistry. Reagent strip is almost made from cellulose. Reagent areas are coated with chemicals which react with constituents in urine. A change of color depends on the amount of the constituents. The result values are demonstrated in Table 1. There are

single and multiple reagent strips for urinalysis. The machine has been widely used in health check-up. Here, the authors reported a comparative study between two automated urine strip analyzer machines from different companies used at King Chulalongkorn Memorial Hospital. Good agreement between both machines can be derived.

Materials and Methods

This study was comparative between two automated urine strip analyzers: Miditron M and Uriscan Pro⁺. ^(5,6) Routine urinalysis was analyzed by Miditron M. YD Diagnostics support Uriscan Pro⁺ to compare with Miditron M. One hundred and ninety-nine individual urine samples were included (Table 2). For each subject, the urine strip tests by both methods were performed.

Table 1. Results from urine strip reader (5,6)

Parameters	Results from urine strip reader			
PH	5 5.5 6 6.5 7 8			
Glucose	Negative 1+ 2+ 3+ 4+			
Specific gravity	1.000 1.005 1.010 1.015 1.020 1.030			
Bilirubin	Negative 1+ 2+ 3+ 4+			

Table 2. Population classified by sex, age.

Age (year)	Female	Male	Total 13
Under 11	6	7	
11-20	9	2	11
21-30	20	8	28
31-40	32	13	45
41-50	21	9	30
51-60	17	8	25
61-70	15	10	25
Over 70	11	11	22
Total	131	68	199

Concerning Miditron M (Roche Diagnostics), the measurement method is reflectance photometer with color image sensor reader head using two main wavelengths at 447 and 620 nm. ⁽⁷⁾ Concerning Uriscan Pro+ (YD Electronics), the measurement method is reflectance photometer with CCD color image sensor reader head using light wavelength visible rays ranging from 440 to 700nm. ⁽⁸⁾ All analyses in this study were performed at room temperature according to the machine manufacturers' guidelines under routine laboratory quality control processes.

The parameters used for comparative study in this study included pH, glucose, specific gravity and bilirubin. The correlations between the biochemical parameter results from both analyzers were assessed using percentage of agreement of the results, and interclass correlation co-efficient.

Results

All 199 samples were analyzed by both Miditron M and Uriscan Pro+. Concerning the correlation between both machines, the results are shown in Table 3, Table 4.

Table 3. Correlation between the biochemical parameter results from Miditron M and Uriscan Pro+.

Agreement (%	
70.8 %	
100 %	
65.3 %	
100 %	

Table 4. Intraclass correlation coefficient 95 % Confidence Interval pH and specific gravity results from Miditron M and Uriscan Pro+.

Parameters	Intraclass correlation coefficient	95% Confidence Interval		
		Lower	Upper	p-value
		Bound	Bound	
PH	0.89	0.85	0.91	<0.001
Specific gravity	0.82	0.76	0.86	<0.001

Discussion

Urine biochemical test is a basic test practiced in laboratory medicine. With advanced automated analyzer technology, many new automated urine strip analyzers have been produced. (4) Based on laboratory medicine principles, it is necessary to evaluate the different machines for analysis of the same test since there might be certain degree of significant discrimination among those machines. Here, the author performed a comparative study on two automated urine strip analyzers. Both analyzers made use of the same reflectance photometer principles with different wavelengths detection.

According to this study, the author found that the two machines did provide 100 % agreement of glucose and bilirubin results from both machines. Although the agreement of pH and specific gravity determined by both machines equals 70.8 % and 65.3 %, respectively, the correlation coefficients between pH, specific gravity using Miditron M and Uriscan Pro⁺ were 0.89 and 0.82; p < 0.001. The authors hereby conclude that the use of either urine strip analyzer machine can provide similar analytical results for only glucose and bilirubin. The interpretation of urine pH and specific gravity should be careful. Confirmation by other standard methods such as spectrometry for urine specific gravity is needed.

Acknowledgements

We would like to thank Dr. Junya Pattaraarchachai of the Graduate Studies Faculty of Medicine, Thammasat University for her help in data analysis and every staff member of the Division of Laboratory Medicine, King Chulalongkorn Memorial Hospital for their good cooperation.

References

- Free AH, Free HM. Urinalysis , critical discipline of clinical science. CRC Crit Rev Clin Lab Sci 1972 Dec; 3(4): 481-531
- Brunzel NA. Fundamentals of Urine and body fluid analysis. Pensulvania: W.B. Saunders, 1994.
- Wiwanitkit V. Abnormal laboratory results as presentation in screening test. Chula Med J 1998 Dec; 42(12): 1059-67
- 4.King C. Automated methods in urinalysis. Clin Lab Sci 1998 Jan-Feb; 11(1): 44-5
- Roche Diagnostics Ltd. Coubur 10 test M, reagent strip manual. Bell Lane, BN71LG, Lewes, E. Sussex, GB
- YD Diagnostics. Uriscan. Reagent strips manual. Seoul, Korea
- Roche Diagnostics Ltd. Miditron Operator's Manual. Bell Lane, BN7 1LG, Lewes, E. Sussex, GB
- 8. YD Electronic Co., Ltd. Uriscan Pro Urineanalyzer
 Operating Manual. Seoul ,Korea